

*Methods in Enzymology*

*Volume 404*

# *GTPases Regulating Membrane Dynamics*

EDITED BY

*William E. Balch*

DEPARTMENT OF CELL BIOLOGY  
THE SCRIPPS RESEARCH INSTITUTE  
LA JOLLA, CALIFORNIA

*Channing J. Der*

DEPARTMENT OF PHARMACOLOGY  
THE UNIVERSITY OF NORTH CAROLINA AT CHAPEL HILL  
CHAPEL HILL, NORTH CAROLINA

*Alan Hall*

CRC ONCOGENE AND SIGNAL TRANSDUCTION GROUP  
MRC LABORATORY FOR MOLECULAR CELL BIOLOGY  
UNIVERSITY COLLEGE LONDON  
LONDON, UNITED KINGDOM

Universitäts- und Landes-  
bibliothek Darmstadt  
Bibliothek Biologie

Div. für. .... 15.859  
.....



AMSTERDAM • BOSTON • HEIDELBERG • LONDON  
NEW YORK • OXFORD • PARIS • SAN DIEGO  
SAN FRANCISCO • SINGAPORE • SYDNEY • TOKYO

Academic Press is an imprint of Elsevier



## Table of Contents

CONTRIBUTORS TO VOLUME 404 . . . . .	xi
PREFACE . . . . .	xix
VOLUMES IN SERIES . . . . .	xxi
1. Arrays of Transfected Mammalian Cells for High Content Screening Microscopy	HOLGER ERFLE AND RAINER PEPPERKOK 1
2. Imaging Platforms for Measurement of Membrane Trafficking	RAINER PEPPERKOK, JEREMY C. SIMPSON, JENS RIETDORF, CIHAN CETIN, URBAN LIEBEL, STEFAN TERJUNG, AND TIMO ZIMMERMANN 8
3. Bioinformatic “Harvester”: A Search Engine for Genome-Wide Human, Mouse, and Rat Protein Resources	URBAN LIEBEL, BJOERN KINDLER, AND RAINER PEPPERKOK 19
4. Microinjection as a Tool to Explore Small GTPase Function	BRIAN STORRIE 26
5. Visualizing Intracellular Events <i>In Vivo</i> by Combined Video Fluorescence and 3-D Electron Microscopy	ALEXANDER A. MIRONOV, GALINA V. BEZNOUSSENKO, ALBERTO LUINI, AND ROMAN S. POLISHCHUK 43
6. Recombinant Production in Baculovirus-Infected Insect Cells and Purification of the Mammalian Sec13/Sec31 Complex	CEMAL GURKAN AND WILLIAM E. BALCH 58
7. Purification and Properties of Mammalian Sec23/24 from Insect Cells	PAUL LAPOINTE AND WILLIAM E. BALCH 66
8. Purification and Functional Properties of Yeast Sec12 GEF	EUGENE FUTAI AND RANDY SCHEKMAN 74
9. Reconstitution of Cargo-Dependent COPII Coat Assembly on Proteoliposomes	KEN SATO AND AKIHICO NAKANO 83
10. Real-Time Assays for the Assembly-Disassembly Cycle of COP Coats on Liposomes of Defined Size	JOËLLE BIGAY AND BRUNO ANTONNY 95

11. Assay and Measurement of Phospholipase D Activation by Sar1	MEIR ARIDOR, KUNTALA SHOME, AND GUILLERMO ROMERO	108
12. Use of Acyltransferase Inhibitors to Block Vesicular Traffic Between the ER and Golgi Complex	WILLIAM J. BROWN AND JOHN A. SCHMIDT	115
13. Tethering Assays for COPI Vesicles Mediated by Golgins	AYANO SATOH, JÖRG MALSAM, AND GRAHAM WARREN	125
14. Investigating the Role of ADP-Ribosylation Factor 6 in Tumor Cell Invasion and Extracellular Signal-Regulated Kinase Activation	HOLLY HOOVER, VANDHANA MURALIDHARAN-CHARI, SARAH TAGUE, AND CRISLYN D'SOUZA-SCHOREY	134
15. Assays and Properties of the Arf GAPs AGAP1, ASAP1, and Arf GAP1	MAGNUS MUTAH CHE, ZHONGZHEN NIE, AND PAUL A. RANDAZZO	147
16. Preparation of Myristoylated Arf1 and Arf6	VI LUAN HA, GERAINT M. H. THOMAS, STACEY STAUFFER, AND PAUL A. RANDAZZO	164
17. BIG1 and BIG2, Brefeldin A-Inhibited Guanine Nucleotide-Exchange Factors for ADP-Ribosylation Factors	HEATHER D. JONES, JOEL MOSS, AND MARTHA VAUGHAN	174
18. Cytohesin-1: Structure, Function, and ARF Activation	GUSTAVO PACHECO-RODRIGUEZ, JOEL MOSS, AND MARTHA VAUGHAN	184
19. ADP-Ribosylation Factor Domain Protein 1 (ARD1), a Multifunctional Protein with Ubiquitin E3 Ligase, GAP, and ARF Domains	ALESSANDRO VICHI, JOEL MOSS, AND MARTHA VAUGHAN	195
20. Expression of BIG2 and Analysis of Its Function in Mammalian Cells	HYE-WON SHIN, CHISA SHINOTSUKA, AND KAZUHISA NAKAYAMA	206
21. Assays and Properties of the ArfGAPs, AMAP1 and AMAP2, in Arf6 Function	SHIGERU HASHIMOTO, ARI HASHIMOTO, ATSUKO YAMADA, YASUHITO ONODERA, AND HISATAKA SABE	216

---

22. Functional Assay of EFA6A, a Guanine Nucleotide Exchange Factor for ADP-Ribosylation Factor (ARF6), in Dendritic Formation of Hippocampal Neurons	HIROYUKI SAKAGAMI, AKIFUMI KAMATA, KOHJI FUKUNAGA, AND HISATAKE KONDO	232
23. Functional Assay of ARNO and ARF6 in Neurite Elongation and Branching	DELIA J. HERNÁNDEZ-DEVIEZ AND JEAN M. WILSON	242
24. Analysis of the Interaction Between Cytohesin 2 and IPCEF1	KANAMARLAPUDI VENKATESWARLU	252
25. Assay and Properties of the GIT1/p95-APP1 ARFGAP	IVAN DE CURTIS AND SIMONA PARIS	267
26. Preparation and Characterization of Recombinant Golgin Tethers	AYANO SATOH, MATTHEW BEARD, AND GRAHAM WARREN	279
27. Purification and Functional Properties of the Membrane Fissioning Protein CtBP3/BARS	CARMEN VALENTE, STEFANIA SPANÒ, ALBERTO LUINI, AND DANIELA CORDA	296
28. <i>In Vitro</i> Assays of Arf1 Interaction with GGA Proteins	HYE-YOUNG YOON, JUAN S. BONIFACINO, AND PAUL A. RANDAZZO	316
29. The Role of EFA6, Exchange Factor for Arf6, for Tight Junction Assembly, Functions, and Interaction with the Actin Cytoskeleton	FRÉDÉRIC LUTON	332
30. <i>In Vitro</i> Reconstitution of ARF-Regulated Cytoskeletal Dynamics on Golgi Membranes	JI-LONG CHEN, WEIDONG XU, AND MARK STAMNES	345
31. Assays and Properties of Arfaptin 2 Binding to Rac1 and ADP-Ribosylation Factors (Arfs)	OK-HO SHIN AND JOHN H. EXTON	359
32. Analysis of Arf Interaction with GGAs <i>In Vitro</i> and <i>In Vivo</i>	KAZUHISA NAKAYAMA AND HIROYUKI TAKATSU	367
33. Arf6 Modulates the $\beta$ -Actin Specific Capping Protein, $\beta$ cap73	ALICE Y. WELCH, KATHLEEN N. RILEY, CRISLYN D'SOUZA-SCHOREY, AND IRA M. HERMAN	377
34. Functional Assay of Effectors of ADP Ribosylation Factor 6 During Clathrin/AP-2 Coat Recruitment to Membranes	MICHAEL KRAUSS AND VOLKER HAUCKE	388

35. Assays to Study Phospholipase D Regulation by Inositol Phospholipids and ADP-Ribosylation Factor 6	DALE J. POWNER, TREVOR R. PETTITT, AND MICHAEL J. O. WAKELAM	398
36. Assay and Functional Properties of the Tyrosine Kinase Pyk2 in Regulation of Arf1 Through ASAP1 Phosphorylation	ANAMARIA KRULJAC-LETUNIC AND ANDREE BLAUKAT	411
37. ADP-Ribosylation Factor 6 Regulation of Phosphatidylinositol-4,5-Bisphosphate Synthesis, Endocytosis, and Exocytosis	YOSHIKATSU AIKAWA AND THOMAS F. J. MARTIN	422
38. Interaction of Arl1 GTPase with the GRIP Domain of Golgin-245 as Assessed by GST (Glutathione-S-transferase) Pull-Down Experiments	LEI LU, GUIHUA TAI, AND WANJIN HONG	432
39. Functional Analysis of Arl1 and Golgin-97 in Endosome-to-TGN Transport Using Recombinant Shiga Toxin B Fragment	GUIHUA TAI, LEI LU, LUDGER JOHANNES, AND WANJIN HONG	442
40. Assays Used in the Analysis of Arl2 and Its Binding Partners	J. BRADFORD BOWZARD, J. DANIEL SHARER, AND RICHARD A. KAHN	453
41. Assay and Functional Analysis of the ARL3 Effector RP2 Involved in X-Linked Retinitis Pigmentosa	R. JANE EVANS, J. PAUL CHAPPLE, CELENE GRAYSON, ALISON J. HARDCastle, AND MICHAEL E. CHEETHAM	468
42. Reconstitution of Transport to Recycling Endosomes <i>In Vitro</i>	RENÉ BARTZ, CORINNE BENZING, AND OLIVER ULLRICH	480
43. Robust Colorimetric Assays for Dynamin's Basal and Stimulated GTPase Activities	MARILYN LEONARD, BYEONG DOO SONG, RAJESH RAMACHANDRAN, AND SANDRA L. SCHMID	490
44. Clathrin-Coated Vesicle Formation from Isolated Plasma Membranes	ISHIDO MIWAKO AND SANDRA L. SCHMID	503
45. Nucleotide Binding and Self-Stimulated GTPase Activity of Human Guanylate-Binding Protein 1 (hGBP1)	SIMONE KUNZELMANN, GERRIT J. K. PRAEFCKE, AND CHRISTIAN HERRMANN	512
46. Stimulation of Dynamin GTPase Activity by Amphiphysin	YUMI YOSHIDA AND KOHJI TAKEI	528

---

47. Tuba, A GEF for CDC42, Links Dynamin to Actin Regulatory Proteins	GIANLUCA CESTRA, ADAM KWIAKOWSKI, MARCO SALAZAR, FRANK GERTLER, AND PIETRO DE CAMILLI	537
48. Expression and Properties of Sorting Nexin 9 in Dynamin-Mediated Endocytosis	RICHARD LUNDMARK AND SVEN R. CARLSSON	545
49. Rapid Purification of Native Dynamin I and Colorimetric GTPase Assay	ANNIE QUAN AND PHILLIP J. ROBINSON	556
50. Assays and Functional Properties of Auxilin-Dynamin Interactions	SANJA SEVER, JESSE SKOCH, BRIAN J. BACSKAI, AND SHERRI L. NEWMYER	570
51. Assay and Functional Analysis of Dynamin-Like Protein 1 in Peroxisome Division	MICHAEL SCHRADER AND STEPHEN J. GOULD	586
52. <i>In Vitro</i> Reconstitution of Discrete Stages of Dynamin-Dependent Endocytosis	MATTHEW K. HIGGINS AND HARVEY T. McMAHON	597
53. A Continuous, Regenerative Coupled GTPase Assay for Dynamin-Related Proteins	ELENA INGERMAN AND JODI NUNNARI	611
54. Assay and Properties of the Mitochondrial Dynamin Related Protein Opa1	LORENA GRIPARIC AND ALEXANDER M. VAN DER BLIEK	620
55. Assay and Functional Analysis of Dynamin-Like Mx Proteins	GEORG KOCHS, MIKE REICHELT, DGANIT DANINO, JENNY E. HINSHAW, AND OTTO HALLER	632
AUTHOR INDEX . . . . .		645
SUBJECT INDEX . . . . .		677